



Clamps and forceps



BACKHAUS

Towel clamp with sharp overlapping tips.

* 107-10910- 11 cm 4.25"



LAHEY

Tissue holding forceps with 3x3 teeth.

* 107-13030- 3x3 Teeth 16 cm 6.25"



107-13030-




ALLIS

Tissue-grasping forceps for secure grip of tissue. Useful in rhytidectomy.

* 107-13120-15 5x6 dull teeth 15 cm 6"
107-13120-21 5x6 dull teeth 21 cm 8.25"



107-13120-15
107-13120-21



C
L
A
M
P
S
&
F
O
R
C
E
P
S


* Photo shows this instrument. Not shown at actual size.



ALLIS NON-TRAUMATIC

Atraumatic tissue-grasping forceps. Useful in rhytidectomy.



* 107-13130- Atraumatic 20 cm 8"



FOERSTER

Sponge holding forceps with several different loops.

- 107-12485- Small, straight Serrated loops 24 cm 9.5"
- 107-12486- Small, curved Serrated loops 24 cm 9.5"
- 107-12487- Large, straight Serrated loops 24 cm 9.5"
- 107-12488- Large, curved Serrated loops 24 cm 9.5"



CLAMPS & FORCEPS

* Photo shows this instrument. Not shown at actual size