




Stille micro instruments

Stille micro instruments

All edges, the spring unit, pivot and the grip are designed so that gloves cannot be damaged.

The handle has been lengthened up to the pivot to create a perfectly balanced instrument permitting a good hold.



Grooved handle in order to reduce the weight of the instrument.

Anatomically adapted to the surgeon's hand.

The satin finished surface guarantees freedom from glare. The surgeon's eyes do not become tired even during long operations.


Stille manufactures a line of ultra high quality micro instruments, designed to facilitate the surgical procedure. For this purpose they all have the same length and weight. In addition they all require the same amount of pressure for smooth opening and closing of the instrument. The surgeon does not have to modify or adjust his grip when changing instruments.



STILLE MICRO SCISSORS

Well balanced and ergonomic micro scissors of ultra high quality. Available with normal, Super Cut or serrated blades.

- 117-A00- Straight 16 cm 6.25"
- 117-A01- Straight Super Cut 16 cm 6.25"
- 117-A05- Straight serrated 16 cm 6.25"
- * 117-A10- Curved 16 cm 6.25"
- 117-A11- Curved Super Cut 16 cm 6.25"




STILLE MICRO FORCEPS

Well balanced and ergonomic micro forceps of ultra high quality. Straight or curved jaws.




- * 117-B00- Straight 16 cm 6.25"
- 117-B05- Straight Delicate 16 cm 6.25"
- 117-B15- Curved Delicate 16 cm 6.25"
- 117-B20- Straight 1x2 Fine teeth 16 cm 6.25"



STILLE MICRO NEEDLE GUIDE FORCEPS

Well balanced and ergonomic micro needle guide forceps of ultra high quality.

- * 117-B01- Straight Fine ring 16 cm 6.25"



MICRO INSTRUMENTS


* Photo shows this instrument. Not shown at actual size




STILLE MICRO NEEDLEHOLDER

Well balanced and ergonomic micro needle holder. Round handles facilitate finger tip control. Available with straight or curved jaws. Without ratchet.

- * 117-C00- Straight 16 cm 6.25"
- 117-C10- Curved 16 cm 6.25"



117-C00-




117-C10-


STILLE-HEDÉN MICRO NEEDLEHOLDER

Well balanced and ergonomic micro needle holder with delicate tips for fine sutures. Round handles facilitate finger tip control. Available with straight or curved jaws. Without ratchet.

- * 117-C05- Straight Delicate 16 cm 6.25"
- 117-C15- Curved Delicate 16 cm 6.25"



117-C05



117-C15-

M
I
C
R
O

I
N
S
T
R
U
M
E
N
T
S

* Photo shows this instrument. Not shown at actual size.



Instrument care & Single-use products



STILLE INSTRUMENT OIL

Oil for surgical instruments, package of 2 pcs á 50 ml.

* 119-78900-1 2 x 50 ml



INSTRUMENT CARE

STILLE ORIGINAL CLEANING BRUSH

Cleaning brush for surgical instruments. Special stainless steel fine threaded, bristle diameter 0.12 mm, height 20 mm. Handy grip of plastic material. Autoclavable at 121C.

* 119-78950- 20 cm 8" Package of 2 pcs



* Photo shows this instrument. Not shown at actual size.



INSTRUMENT STERILIZATION TRAYS

Instrument trays for transport, sterilization and storage of delicate instruments. The trays are manufactured with a unique "pebble-like" interior surface to assure proper penetration of sterilizing agent between the mat and the tray surface. High quality latches to ensure long lifetime. Silicone mat included. Available in three different sizes. The trays are made to withstand common sterilization methods, including steam, gas, cold sterilizing solutions, and temperatures up to 182 Celsius / 360 Fahrenheit.

199-60010-	Micro Tray	10.1x19x1.9 cm	4x7.5x0.75"
199-60020-	Micro Tray	15.2x25.4x1.9 cm	6x10x0.75"
199-60030-	Instrument Tray	25.4x38.1x1.9 cm	10x15x0.75"



I
N
S
T
R
U
M
E
N
T

C
A
R
E

* Photo shows this instrument. Not shown at actual size



TIP PROTECTORS

Protect instruments from costly repairs. Available in a variety of sizes and colours. Round and flat caps. All tips are latex free, radiopaque and steam sterilizable. Also available with air slits to prevent vacuum seal. Delivered as non-sterile.

199-80010-	Tinted White	1.6x19 mm
199-80015-	Tinted Blue	2x19 mm
199-80020-	Tinted Green	2.8x19 mm
199-80025-	Tinted Red	3.2x25.4 mm
199-80030-	Tinted yellow	4.8x25.4 mm
199-80035-	Tinted Purple	10.3x19 mm
199-80045-	Tinted Flat Orange	1.6x9.5x25.4 mm
199-80040-	Tinted Flat Brown	3.2x15.9x25.4 mm
199-80050-	Tinted Flat Black	3.2x25.4x25.4 mm
199-80055-	Clear with Slits	1.6x19 mm
199-80060-	Clear with Slits	2x19 mm
199-80065-	Clear with Slits	2.8x19 mm
199-80070-	Clear with Slits	3.2x25.4 mm
199-80075-	Clear with Slits	4.8x 25.4 mm
199-80080-	Clear with Slits	10.3x19 mm
199-80085-	Clear with Slits Flat	1.6x9.5x25.4 mm
199-80090-	Clear with Slits Flat	1.6x15.9x25.4 mm
199-80095-	Clear with Slits Flat	3.2x25.4x25.4 mm

T
I
P
P
R
O
T
E
C
T
O
R
S



* Photo shows this instrument. Not shown at actual size.



SURGICAL MARKING PENS

Surgical marking pens in 3 different models. Large ink capacity. Point design allows pen to be angled to surgeon's requirements. Each sterile package includes plastic, flexible ruler that will not rip, tear or smear upon contact with liquids. 18-24 month shelf life from date of manufacture. Delivered in packages of 10 markers.

199-70010-	Violette	Universal Tip	Wide Body
199-70020-	Violette	Ultra Fine Tip	Accu-trim Barrel
199-70030-	Violette	Universal Tip	Accu-trim Barrel



MARKING
PENS

* Photo shows this instrument. Not shown at actual size